

Endocavity Complete Probe Reprocessing Infection Prevention Rounding Checklist

OBSERVATION HEADER
Date:
Department:
Survey Team:
Staff Interviewed:
Equipment Observed:
Notes/Instructions:

1. Staff can show where IFUs are accessed and how they are built into the workflow, as needed.
Compliant? <input type="radio"/> YES or <input type="radio"/> NO
2. PPE is readily and easily accessible where reprocessing activities are taking place.
Compliant? <input type="radio"/> YES or <input type="radio"/> NO
3. Hand hygiene products are available near where the work is being performed and easily accessible.
Compliant? <input type="radio"/> YES or <input type="radio"/> NO
POINT OF USE CLEANING
4. Probe is immediately wiped down after the procedure to remove probe cover, gel, blood, body fluid and to prevent biofilm with an enzymatic product or disinfectant.
Compliant? <input type="radio"/> YES or <input type="radio"/> NO
INSPECTION
5. Staff can demonstrate how they would inspect the probe for damage, including cracks, degradation, tears and sharp edges.
Compliant? <input type="radio"/> YES or <input type="radio"/> NO
6. Staff can verbalize the process if they feel, after inspection, the probe needs to be taken out of service.
Compliant? <input type="radio"/> YES or <input type="radio"/> NO
SOILED TRANSPORT
7. Endocavity or surface ultrasound probe is transported in a sealed container or bag that separates the probe (HLD) and the cord & electrical pack (LLD).
Compliant? <input type="radio"/> YES or <input type="radio"/> NO or <input type="radio"/> N/A (choose N/A if cleaning and disinfection occur in same room as procedure.)
***If automating cleaning and HLD skip to Low-Level Disinfection

MANUAL CLEANING

The following steps are taken per the IFUs:

8. Correct temperature is maintained for the duration of soaking.

Compliant? YES or NO

9. Probe is soaked for the correct amount of time.

Compliant? YES or NO

10. Correct amount of detergent to water is observed.

Compliant? YES or NO

11. Probe is positioned in a way that only the section that should be under water is under the water/detergent solution.

Compliant? YES or NO

12. Rinsing is performed per the IFUs.

Compliant? YES or NO

13. Probes are dried with a single-use, single-packaged low-linting cloth after cleaning is complete.

Compliant? YES or NO

14. Staff dons and doffs appropriate PPE during manual cleaning, per hospital policy, correctly and safely.

Compliant? YES or NO

MANUAL HIGH-LEVEL DISINFECTION

The following steps are taken per the IFUs:

15. Correct temperature is maintained for the duration of soaking

Compliant? YES or NO

16. MRC testing is performed via test strip and verified.

Compliant? YES or NO

17. MRC testing result is documented for all probes undergoing HLD.

Compliant? YES or NO

18. Expiration date on test strip bottle is observed and validated.

Compliant? YES or NO

19. Staff dons and doffs appropriate PPE during manual cleaning, per hospital policy, correctly and safely.

Compliant? YES or NO

20. Probe is soaked for the correct amount of time

Compliant? YES or NO

21. Probe is positioned in a way that only section that should be underwater is under the water/disinfectant solution

Compliant? YES or NO

22. Rinsing per the IFU
Compliant? <input type="radio"/> YES or <input type="radio"/> NO
23. Probes are dried with a single-use, single-packaged low-linting cloth after removal from the HLD machine or manual HLD solution.
Compliant? <input type="radio"/> YES or <input type="radio"/> NO
24. Documentation is observed to verify all steps of the manual HLD cycle have been performed and were either successful or failed.
Compliant? <input type="radio"/> YES or <input type="radio"/> NO
25. Traceability is verified for patient identification.
Compliant? <input type="radio"/> YES or <input type="radio"/> NO
LOW-LEVEL DISINFECTION:
26. Transducer cord, handle and electrical pack are low-level disinfected with a hospital-approved disinfectant.
Compliant? <input type="radio"/> YES or <input type="radio"/> NO
27. Staff can speak to the contact time of the disinfectant used in the prior step.
Compliant? <input type="radio"/> YES or <input type="radio"/> NO
AUTOMATED CLEANING/HIGH-LEVEL DISINFECTION
The following steps are taken per the IFUs:
28. Probes are dried with a single-use, single-packaged low-linting cloth after removal from the cleaner/HLD machine.
Compliant? <input type="radio"/> YES or <input type="radio"/> NO
29. Documentation is observed to verify all steps of the automated cleaning/HLD cycle have been performed and either were successful or failed.
Compliant? <input type="radio"/> YES or <input type="radio"/> NO
30. Traceability is verified for patient identification.
Compliant? <input type="radio"/> YES or <input type="radio"/> NO
DRYING
31. Probes are dried with a single-use, single-packaged low-linting cloth after removal from the HLD machine or manual HLD solution.
Compliant? <input type="radio"/> YES or <input type="radio"/> NO
STORAGE
32. Endocavity or surface ultrasound probe is hung vertically without touching anything. Storage should be in a clean, well-ventilated space. (Check the bottom of the cabinet for stains or towels/pads to indicate probes were not dry when stored.)
Compliant? <input type="radio"/> YES or <input type="radio"/> NO
33. Staff can speak to the process and frequency of cleaning the probe storage cabinet.
Compliant? <input type="radio"/> YES or <input type="radio"/> NO
34. Storage cabinet filter change is performed per IFU with documentation present.
Compliant? <input type="radio"/> YES or <input type="radio"/> NO

35. Staff can verbalize the maximum hang time or storage time per facility policy before reprocessing in between disinfection events.
Compliant? <input type="radio"/> YES or <input type="radio"/> NO
ENVIRONMENTAL
36. Work flows from dirty to clean.
Compliant? <input type="radio"/> YES or <input type="radio"/> NO
37. Staff can verbalize what to do in the event of a chemical spill in the reprocessing space (check status of spill kit).
Compliant? <input type="radio"/> YES or <input type="radio"/> NO
38. Ceiling tiles are free from stains and damage.
Compliant? <input type="radio"/> YES or <input type="radio"/> NO
39. Flooring is in good repair, free from cracks and damage.
Compliant? <input type="radio"/> YES or <input type="radio"/> NO
40. Reprocessing space shows evidence of negative pressure airflow.
Compliant? <input type="radio"/> YES or <input type="radio"/> NO
MISCELLANEOUS
41. Staff can verbalize how they would trace a probe back to a patient if needed.
Compliant? <input type="radio"/> YES or <input type="radio"/> NO
42. Competencies are performed annually by a qualified person, as well as when there is new equipment, supplies or processes.
Compliant? <input type="radio"/> YES or <input type="radio"/> NO



2179 East Lyon Station Road • Creedmoor, NC 27522
 Phone: +1 (919) 255-9472 • Toll Free: 877-255-9472 • Fax: +1 (919) 528-3400
 csmedicalllc.com